

WARD'S

Simulated ABO Blood Typing

Student Study and Analysis Sheets

Note: In response to the concerns of many science educators, WARD'S has developed an alternative blood typing activity that does not use real blood. Students will follow the same procedure used to type actual human blood and obtain results that closely approximate those obtained by real blood typing.

Introduction

Around 1900, Karl Landsteiner discovered that there are at least four different kinds of human blood, determined by the presence or absence of specific agglutinogens (agglutinating antigens) on the surface of red blood cells (erythrocytes). These antigens have been designated as A and B. Antibodies against antigens A or B begin to build up in the blood plasma shortly after birth, the levels peak at about eight to 10 years of age, and the antibodies remain, in declining amounts, throughout the rest of a person's life. The stimulus for antibody production is not clear; however, it had been proposed that antibody production is initiated by minute amounts of A and B antigens that may enter the body through food, bacteria, or other means. Humans normally produce antibodies against those antigens that are not on their erythrocytes: A person with A antigens has anti-B antibodies; a person with B antigens has anti-A antibodies; a person with neither A or B antigens has both anti-A and anti-B antibodies; and a person with both A and B antigens has neither anti-A nor anti-B antibodies: Blood type is based on the antigens, not the antibodies, a person possesses.

The four blood groups are types A, B, AB, and O. Blood type O, characterized by the absence of A or B agglutinogens, is the most common in the United States, in 4596 of the population. Type A is next in frequency, found in 39% of the population. The incidences of types B and AB are 1296 and 496 respectively.

ABO System

Blood Type	Antigens on Erythrocytes (Agglutinogens)	Antibodies in Plasma (Agglutinins)	Can Give Blood To	Can Receive Blood From
A	A	Anti-B	A, AB	O, A
B	B	Anti-A	B, AB	O, B
AB	A and B	Neither Anti-A nor Anti-B	AB	O, A, B, AB
O	Neither A nor B	Both Anti-A and Anti-B	O, A, B, AB	O

Process of Agglutination

There is a simple test to determine blood type, performed with antisera containing high levels of anti-A and anti-B agglutinins. Several drops of each kind of antiserum are added to separate samples of blood. If agglutination (clumping) occurs only in the suspension to which the anti-A serum was added, the blood type is A. If agglutination occurs only in the anti-B mixture, the blood type is B. Agglutination in both samples indicates that the blood type is AB. The absence of agglutination in any sample indicates that the blood type is O.

Agglutination Reaction of ABO Blood-Typing Sera

Reaction		
Anti-A Serum	Anti-B Serum	Blood Type
Agglutination	No Agglutination	Type A
No Agglutination	Agglutination	Type B
Agglutination	Agglutination	Type AB
No Agglutination	No Agglutination	Type O

Importance of Blood Typing

As noted in the table above, people can receive transfusions of only certain blood types, depending on the type of blood they have. If incompatible blood types are mixed, erythrocyte destruction, agglutination and other problems can occur. For instance, if a person with Type B blood is transfused with blood type A, the recipient's anti-A antibodies will attack the incoming type A erythrocytes. The type A erythrocytes will be agglutinated, and hemoglobin will be released into the plasma. In addition, incoming anti-B antibodies of the type A blood may also attack the type B erythrocytes of the recipient, with similar results. This problem may not be serious, unless a large amount of blood is transfused.

The ABO blood groups and other inherited antigen characteristics of red blood cells are often used in medicolegal situations involving identification of disputed paternity. A comparison of the blood groups of mother, child, and alleged father may exclude the man as a possible parent. Blood typing does not prove that an individual is the father of a child; it merely indicates whether or not he is a possible parent. For example, a child with a blood type of AB, whose mother is type A, could not have as a father a man whose blood type is O.

The Genetics of Blood Types

The human blood types (A, B, AB, and O) are inherited by multiple alleles—three or more genes that occupy a single locus on a chromosome. Gene I^A codes for the synthesis of antigen (agglutinogen) A, gene I^B codes for the production of antigen B on the red blood cells, and gene i (i^o) does not produce any antigens. The phenotypes listed in the table below are produced by the combinations of the three different alleles: I^A , I^B , i^o . When genes I^B and I^A are present in an individual, both are fully expressed. Both I^A and I^B are dominant over i^o ; the genotype of an individual with blood type O must be $i^o i^o$.

Phenotypes	Possible Genotypes
A	$I^A I^A$ $I^A i$ (or i^o)
B	$I^B I^B$ $I^B i$ (i^o)
O	ii ($i^o i^o$)

Use I^A for antigen A, I^B for antigen B, i or i^o for no antigens present. Genes I^A and I^B are dominant over i (i^o).

AB blood type results when both genes I^A and I^B are present.

Artificial Blood

At times it is difficult to find a correct match for a blood type of a person requiring a transfusion. It would be ideal to have some type of artificial blood or blood substitute that wouldn't need to be matched to a patient's blood type; it could save thousands of lives each year. Although the research for artificial blood and blood substitutes continues, it may take years before one is available.

In 1966, Dr. Leland C. Clark, of the University of Cincinnati's College of Medicine, developed the first artificial blood prototype. This milky white solution, which can carry twice as much oxygen as blood does, is a fluoro-carbon emulsion called Fluosol. It is made up of two fluorocarbons, a number of salts, water, and fine particles that are $\frac{1}{70}$ the size of erythrocytes. Because these particles are so small, they can pass through narrowed arteries that the larger erythrocytes can't get through, making it an ideal blood substitute for heart attack and stroke victims; they would recover faster and have less tissue damage. Fluosol has been approved for use in Canada, Holland, and Italy. The US Food and Drug Administration is reviewing Fluosol for use in the United States.

Anthony Hunt and colleagues at the University of California at San Francisco are working with artificial red blood cells called neohemocytes. Neohemocytes, which are microscopic spheres of hemoglobin surrounded by lipids, are capable of carrying oxygen. These microspheres are proving to be a successful substitute for erythrocytes. As with Fluosol, their small size allows them to pass through restricted vessels that might not allow the passage of erythrocytes.

Objective

To use WARD'S Simulated Blood to determine the blood type of four individuals. Also to use a simplified counting technique to estimate the number of red and white blood cells per cubic millimeter.

A. ABO Blood Typing

Materials Needed per Lab Group

- 4 Blood Typing Slides
- 8 Toothpicks

Shared Materials

- 4 Unknown Blood Samples
 - Mr. Smith
 - Mr. Jones
 - Mr. Green
 - Ms. Brown
- Simulated Anti-A Typing Serum
- Simulated Anti-B Typing Serum

Procedure

1. Label each blood typing slide:
 - Slide #1: Mr. Smith
 - Slide #2: Mr. Jones
 - Slide #3: Mr. Green
 - Slide #4: Ms. Brown
2. Place three to four drops of Mr. Smith's blood in each of the A and B wells of Slide #1.
3. Place three to four drops of Mr. Jones's blood in each of the A and B wells of Slide #2.
4. Place three to four drops of Mr. Green's blood in each of the A and B wells of Slide #3.
5. Place three to four drops of Ms. Brown's blood in each of the A and B wells of Slide #4.
6. Add three to four drops of the simulated anti-A serum in each A well on the four slides.
7. Add three to four drops of the simulated anti-B serum in each B well on the four slides.
8. Stir each mixture with a different clean toothpick. Use only one toothpick per well to avoid cross contamination.
9. Examine each well for agglutination. Agglutination indicates a positive test result. Record results in Table 1.

Analysis

Table 1 Agglutination Reactions

	Anti-A Serum	Anti-B Serum	Blood Type
Slide #1: Mr. Smith			
Slide #2: Mr. Jones			
Slide #3: Mr. Green			
Slide #4: Ms. Brown			

B. (Optional) Blood Cell Count

Note: WARD'S Simulated Blood contains components that simulate red and white blood cells; they are **similar in proportion** and size to those found in real human blood and can be seen under the microscope without staining.

Materials Needed per Lab Group

- 1 Microscope Slide
- 1 Coverslip
- WARD'S Simulated Blood
- Compound Microscope (400X magnification)

Procedure

1. Thoroughly shake one of the vials of WARD'S Simulated Blood. Add a small drop of simulated blood to a microscope slide. Cover with a coverslip.

Note: Lower the coverslip slowly to avoid trapping air bubbles on the slide.

2. Use a paper towel to absorb any excess sample from the edge of the coverslip.
3. Examine the slide with the low power objective. Find the area of the slide with the best distribution of cells-one that has a minimum of 20 red blood cells.
4. Switch to 400X magnification and count the number of simulated red blood cells. Record the number in

Table 2.

Note: Like real blood cells stored at room temperature for long periods of time, some cells will clump together. Count the cells in any clump separately.

5. Count the number of simulated white blood cells. Record the number in Table 2.
6. Repeat the counting procedure with two other fields of view. Record these counts in Table 2.

Analysis

1. Calculate the average of the three cell counts.
2. Multiply the average number of cells by the dilution factor to determine the number of red and white blood cells per cubic millimeter.

Note: The average count is multiplied by the factor indicated in the table for simulated red and white blood cells to correct for dilution and for the fact that only a small volume of blood was observed. The method used is not the actual blood cell counting procedure. A hemacytometer may be used to demonstrate to students the actual counting technique used in clinical settings.

Blood Cell Type	Cell Count			Table 2		Total No. Blood Cells per mm ³ (Avg. No. of Cells x Dilution Factor)
	1	2	3	Avg. No. of Cells	Dilution Factor	
Red (Red)					150,000	
White (Blue)					5,000	

Questions

ABO Blood Group

1. What is Mr. Smith's blood type? What ABO agglutinogens are present in his red blood cells?

2. What is Mr. Green's blood type? What ABO agglutinins are present in the plasma of his blood?

3. What is Mr. Jones's blood type? If Mr. Jones needed a transfusion, what blood type(s) could he safely receive?

9. Could a man with an AB blood type be the father of a child with type O **blood**?

10. Could a man with an O blood type be the father of a child with type AB blood?

11. Could a type B child with a type A mother have a type A father?

12. What are the possible genetic combinations of a child whose parents' blood types are A and B?

Science in Motion

Materials List

Lab: Simulated ABO Blood Typing

Number of Lab Groups Prepared: _____

Equipment per lab group	Delivered	Returned
4 Blood Typing Slides		
Compound Microscope (400x magnification)		

Consumables	Delivered
8 Toothpicks	
Simulated Blood	
Serum for Blood type A (1 bottle)	
Serum for Blood type B (1 bottle)	
Serum for RH factor (1 bottle)	